

DIAGNOSTIC CONSUMABLES

Diagnostic Reagent Sample Kits

As diagnostic tests move closer to the point-of-care, reagent storage and delivery methods become a challenge. J-Pac Medical's lab-on-chip reagent blisters allow reagents to be stored directly on the microfluidic cartridge, enabling the smaller footprint required by point-of-care diagnostic tests.

Contact customer service at 1-603-692-9955 to order sample kits.

About J-Pac Medical

J-Pac Medical is the trusted manufacturing and packaging outsourcing partner to medical device companies seeking to deliver superior quality, improve time-to-market, and simplify the supply chain for single-use medical devices.

| Dim. (Inches) | Sample Blister Fluid Volume | | | |
|---------------------------|-----------------------------|-------------|-------------|---------------|
| | 30 μ l | 130 μ l | 500 μ l | 1,000 μ l |
| Height | 0.108 | 0.175 | 0.26 | 0.268 |
| Diameter | 0.394 | 0.585 | 0.875 | 0.988 |
| Flange | 0.065 | 0.118 | 0.096 | 0.118 |
| Pierce Dia. | 0.264 | 0.375 | 0.625 | 0.752 |
| Quantity | Water Filled (In-Stock) | | | |
| 25 | \$250 | | | |
| 100 | \$500 | | | |
| 250 | \$1,000 | | | |
| 500 | \$1,500 | | | |
| Customer Supplied Reagent | | | | |
| <1000 | \$1500 plus \$5/ea | | | |
| 1000-2000 | \$1500 plus \$4/ea | | | |

