

DIAGNOSTIC CONSUMABLES

Diagnostic Reagent Sample Kits

As diagnostic tests move closer to the point-of-care, reagent storage and delivery methods become a challenge. J-Pac Medical's lab-on-chip reagent blisters allow reagents to be stored directly on the microfluidic cartridge, enabling the smaller footprint required by point-of-care diagnostic tests.

Contact customer service at 1-603-692-9955 to order sample kits. Credit Cards Accepted.

Dimensions	Sample Blister Fluid Volume					
	30 μ l	130 μ l	250 μ l	500 μ l	750 μ l	1,000 μ l
Height	0.108	0.175	0.190	0.26	0.245	0.268
Diameter	0.394	0.585	0.765	0.875	1.006	0.988
Flange	0.065	0.118	0.133	0.096	0.125	0.118
Pierce Dia.	0.264	0.375	0.500	0.625	0.750	0.752
Quantity	Water Filled (In-Stock)					
25	\$250					
100	\$500					
250	\$1,000					
500	\$1,500					
	Customer Supplied Reagent					
<1000	\$1500 plus \$5/ea					
1000-2000	\$1500 plus \$4/ea					

About J-Pac Medical

J-Pac Medical is the trusted manufacturing and packaging outsourcing partner to medical device companies seeking to deliver superior quality, improve time-to-market, and simplify the supply chain for single-use medical devices.

